

GE Healthcare

# Giraffe Warmers

For Neonatal Intensive Care

Helping make sick babies well...



# We want you to notice what we're missing

The Giraffe® Warmer from GE is unlike any warmer you have seen...outstanding thermal performance WITHOUT a clumsy overhead heater!

Performing procedures such as X-rays and surgery just got easier. Attending clinicians remain cool and comfortable. Your baby always stays warm and content.

NICU babies in any environment need a warm, comfortable bed to support their development. The Giraffe Warmer provides uninterrupted thermal support while at the same time giving the caregiver complete access and visibility. Newly integrated features and options minimize the amount of extra equipment in the care area.

The new Hourglass Heat Profile optimizes heat coverage for the mattress and baby while minimizing heat to caregivers. The recessed heater provides exceptional access and consistent heat therapy during interventions.




# Reshaping neonatal care

The controls have a full-colour display that is easy to read, even from a distance. Full colour trends track temperature, pulse rate, SpO<sub>2</sub> and weight. The Hands Free Alarm Silence allows you to control alarms without touching a surface during a procedure. The elimination of alarms during warm-up mode keeps the area quiet and friendly for the baby, family and staff.

The built-in dimmable observation lights and aimable procedure light are standard on the Giraffe Warmer and reduce the need for additional lighting. They are fully adjustable to provide the lighting you need, when and where you need it.

Advanced care options are available for the NICU too. The integrated resuscitation option helps to standardize the resuscitation practice within the NICU. Capture oxygen saturation data with the integrated SpO<sub>2</sub> option. Both options are fully integrated and always ready when you need them.



## Heat where you need it

The Hourglass Heat Profile provides uniform coverage across the mattress.

# Transforming care delivery

## For the clinician...

- Total access
- Complete control
- A New level of neonatal care

## For the baby...

- A warm, comfortable environment
- Nurturing care
- Helping make sick babies well

© 2009 General Electric Company – All rights reserved.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Representative for the most current information.

Giraffe is a trademark of Datex-Ohmeda, Inc.

GE and GE Monogram are trademarks of General Electric Company.

All other company and product names mentioned may be trademarks of the companies with which they are associated.

GE Healthcare Finland Oy., doing business as GE Healthcare.

## About GE Healthcare

GE Healthcare provides transformational medical technologies and services that are shaping a new age of patient care. Our broad expertise in medical imaging and information technologies, medical diagnostics, patient monitoring systems, drug discovery, biopharmaceutical manufacturing technologies, performance improvement and performance solutions services helps our customers to deliver better care to more people around the world at a lower cost. In addition, we partner with healthcare leaders, striving to leverage the global policy change necessary to implement a successful shift to sustainable healthcare systems.

Our “healthymagination” vision for the future invites the world to join us on our journey as we continuously develop innovations focused on reducing costs, increasing access and improving quality and efficiency around the world. For more information about GE Healthcare, visit our website at [www.gehealthcare.com](http://www.gehealthcare.com)

GE Healthcare

P.O. Box 900, FIN-00031 GE, Finland

Tel. +358 10 394 11 • Fax +358 9 146 3310

[www.gehealthcare.com](http://www.gehealthcare.com)



GE imagination at work