

GE Healthcare

Aespire* View

Clarity of view...clarity of decision



Advanced anesthesia delivery in a compact system

You told us you have many different types of patients undergoing different procedures ranging from routine surgery to those with complex trauma. That you need a compact, integrated anesthesia delivery system that helps you deliver attentive care. Finally, you said you need a well-designed, user-friendly system that lets you focus on your patients, not on the equipment.

Working closely with you and listening to what really impacts your practice, GE Healthcare developed the Aespire View anesthesia delivery system.

A cost effective solution for hospital operating rooms and standalone surgery centers, the Aespire View demonstrates how GE's advanced technology leadership is helping drive the future of perioperative care.

Designed for the way you work

Clarity of view

Large 12-inch color ventilator display with color waveform and alarm message indicators along with familiar user interface helps you precisely assess and quickly respond to changes in your patient's status.

Helps manage costs

Total fresh gas or total flow sensing (TFS) module (optional) electronically measures, displays, and exports both individual gas flow and total fresh gas flow delivered by the flow tubes. The additional precision afforded by the electronic measurement helps in managing costs.

	MV l/min	RR /min
	6.1 10.5 0.5	16
Circle	Fresh Gas Flow l/min	
More Settings	O2 3.2	6.1 Total
	Air 2.9	
	N2O 0.00	

Up and running at a moment's notice

The Aespire View fast and easy power up helps you waste less time waiting for your equipment.

Precision flow technology you can count on

Manage your patient's gas delivery with the Aespire View's consistent pneumatic gas flow mixing, individual flow controls, and dual-flow tubes. A minimum O₂ flow of 50 mL/min helps ensure accuracy in low-flow situations, and an auxiliary common gas outlet facilitates low-pressure leak testing.

Keep your OR organized and clutter-free

We designed the Aespire View to be lightweight, compact, and easy to move around. Aespire's sleek, contemporary design is complemented by smooth, durable, easy-to-clean surfaces. You'll appreciate some of the features GE has engineered into the system:

- Advanced cable management arm (optional)
- Extra large work surface area
- Bi-level illumination
- Two roomy storage drawers with optional removable inserts for small items.

The Aespire View offers an integrated standard AGCO and optional integrated features including auxiliary O₂, suction control, CO₂ bypass and multiple monitor integration kits.



More ventilation choices means more tailored care

With five ventilation packages, including GE's innovative Pressure Controlled Ventilation-Volume Guaranteed (PCV-VG) technology, the Aespire View lets you tailor ventilation for each patient.

Ventilation modes are available à la Carte, so you can tailor the clinical capabilities of your anesthesia delivery system with the options you need.



Volume Control Ventilation delivers the tidal volume you set, regardless of changes in the fresh gas flow and airway pressure, up to the selected pressure limit.

Pressure Control Ventilation delivers the maximum tidal volume for the desired time and pressure selected.

Synchronized Intermittent Mandatory Ventilation (SIMV) and Pressure Supported Ventilation (PSVPro*) simplify the care of your spontaneously breathing patients. Both modes offer an adjustable flow trigger, trigger window, electronic PEEP, and an apnea backup is available with PSVPro.

Pressure Controlled Ventilation – Volume Guaranteed (PCV-VG) gives you the best of both worlds: the benefits of pressure controlled ventilation and a consistent tidal volume.

Precise gas delivery

GE has built its Advanced Breathing System (ABS) into the Aespire View. By providing a system well-suited for low and minimal flow anesthesia, the ABS may help reduce patients' agent consumption. Integrated electronics detect when the ABS is completely engaged or disengaged.

- With its small breathing system volume, the ABS responds to fresh gas changes rapidly even during low flow anesthesia
- Fewer parts to help prevent unwanted disconnects and leaks
- Bag-to-vent switch works without electronics – in case of power failure
- Autoclavable and latex-free surfaces
- Optional EZChange and condenser
- Electronic detection of “absorber off” and optional CO₂ bypass



Exceptional support sets us apart

The Aespire View helps keep affordability in reach with:

- Only one scheduled maintenance check per year
- Tec* 7 vaporizers are planned factory service free
- Support for minimal and low-flow anesthesia

GE Healthcare is deeply committed to supporting you long after the point of purchase. Our world class service teams offer you flexible solutions to help ensure that your equipment performs at exceptional standards throughout its lifecycle. Together, we're focused on helping to making sure your investment will serve your long-term requirements.

© 2011 General Electric Company – All Rights Reserved.
General Electric Company reserves the right to make changes in specification and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE representative for the most current information.

GE and GE Monogram are trademarks of General Electric Company.

*Aespire, Tec and PSVPro are trademarks of General Electric Company.

Datex-Ohmeda, Inc., a General Electric company, doing business as GE Healthcare.

About GE Healthcare

GE Healthcare provides transformational medical technologies and services that are shaping a new age of patient care. Our broad expertise in medical imaging and information technologies, medical diagnostics, patient monitoring systems, drug discovery, biopharmaceutical manufacturing technologies, performance improvement and performance solutions services help our customers to deliver better care to more people around the world at a lower cost. In addition, we partner with healthcare leaders, striving to leverage the global policy change necessary to implement a successful shift to sustainable healthcare systems.

Our “healthymagination” vision for the future invites the world to join us on our journey as we continuously develop innovations focused on reducing costs, increasing access and improving quality around the world.

Headquartered in the United Kingdom, GE Healthcare is a unit of General Electric Company (NYSE: GE). Worldwide, GE Healthcare employees are committed to serving healthcare professionals and their patients in more than 100 countries. For more information about GE Healthcare, visit our website at www.gehealthcare.com.

GE Healthcare
3030 Ohmeda Drive
P.O. Box 7550
Madison, WI 53707-7550
U.S.A.

www.gehealthcare.com



imagination at work